




Patient comfort is a priority

Sampling Line Selection Guide

Comfort for your patient, Relaxation for you.

For Intubated and Non-intubated patient



	Endo	Nasal	Nasal O2
Description	Sampling line for COMDEK capnography module		
Order Number	ET-305T-C	NC-005A	NC-004A
Spec.	Adult/Pediatric; standard: 2.5m(8.2ft)	Adult/Pediatric; standard: 3m(9.8ft)	Adult/Pediatric; standard: 3m(9.8ft)
Duration	8hrs	8hrs	8hrs
Packaging	25pcs/ box	25pcs/ box	25pcs/ box
CE	*	*	*
TFDA	*	*	*

For Intubated and Non-intubated patient long-term use



	Endo	Nasal	Nasal O2
Description	Sampling line for COMDEK capnography module		
Order Number	ET-306T-C	NC-006A	NC-007A
Spec.	Adult/Pediatric; standard: 2.5m(8.2ft)	Adult/Pediatric; standard: 3m(9.8ft)	Adult/Pediatric; standard: 3m(9.8ft)
Duration	12hrs	12hrs	12hrs
Packaging	25pcs/ box	25pcs/ box	25pcs/ box
CE	*	*	*
TFDA	*	*	*

For Philips Microstream

For Intubated and Non-intubated patient short-term use



	Endo	Nasal	Nasal O2
Description	Sampling line for PHILIPS Microstream capnography module		
Order Number	ET-301T-P	NC-101A-P	NC-102A-P
Spec.	Adult/Pediatric; standard: 2.5m(8.2ft)	Adult/Pediatric; standard: 3m(9.8ft)	Adult/Pediatric; standard: 3m(9.8ft)
Duration	8hrs	8hrs	8hrs
Packaging	25pcs/ box	25pcs/ box	25pcs/ box
CE	*	*	*
TFDA	*	*	*

For Water Trap

For Intubated and Non-intubated patient long-term use



	Endo	Nasal	Nasal O2
Description	Sampling line for Water Trap capnography module		
Order Number	ST-004	NC-001A	NC-002A
Spec.	Adult/Pediatric; standard: 2.5m(8.2ft)	Adult/Pediatric; standard: 3m(9.8ft)	Adult/Pediatric; standard: 3m(9.8ft)
Packaging	25pcs/ box	25pcs/ box	25pcs/ box
CE	*	*	*
TFDA	*	*	*



COMDEK INDUSTRIAL CORP.

7F., No.1, Alley 6, Lane 235, Pao Chiao Rd.,
Hsin Tien Dist, New Taipei City, 23145, Taiwan

TEL: +886- 2-2918-1266

FAX: +886-2-2918-1277

E-MAIL: info@comdek.com

www.comdek.com



INFO

2025/02/06 V5